Discoveries in IMAGES

Nanoscale imaging of E. coli cells by expansion microscopy

Sharey Cheng, Yongxin Zhao*

Department of Biological Sciences, Carnegie Mellon University, Pittsburgh, PA, USA

**Corresponding authors*: Yongxin Zhao, PhD, Department of Biological Sciences, Carnegie Mellon University, Pittsburgh, PA, USA, Phone: +1 412-268-8236; Email:yongxinz@andrew.cmu.edu.

Submitted: Sept. 27th, 2019; Revised: Sept. 29th, 2019; Accepted: Sept. 29th, 2019; Published: Sept. 30th, 2019; Citation: Cheng S, Zhao Y. Nanoscale imaging of E. coli cells by expansion microscopy. Discoveries 2019, Jul-Sep; 7(3): e98. DOI: 10.15190/d.2019.11

ABSTRACT

Expansion microscopy (ExM) is an emerging superresolution imaging technology. ExM works by biological specimen infusing а with а superabsorbent hydrogel, followed by mechanical homogenization and isotropical expansion of the specimen in water. The unique and cost-effective process of ExM enables super-resolution optical imaging of sample of interest without the need to invest and use of a sophisticated microscope instrument. Here, we demonstrate that a nearly 3fold isotropic physical expansion of E.coli fixed cells can be achieved in PBS, and the cell morphology during binary fission is clearly resolved in the expanded state, using a diffractionlimited microscope.

Keywords

Expansion Microscopy, Escherichia coli, bacteria, super-resolution imaging.

Abbreviations

Expansion microscopy (ExM); Structured illumination microscopy (SIM); Stochastic optical reconstruction Microscopy (STORM); Stimulated emission depletion microscopy (STED); protein retention ExM (proExM); Escherichia coli (E. coli); paraformaldehyde (PFA); DAPI (4',6-diamidino-2phenylindole); wheat germ agglutinin (WGA); numerical aperture (NA); Expansion microscopy (ExM) is a newly emerged imaging technique that seeks to overcome the everlasting hindrance of the optical imaging system, the diffraction limit. Different from the other superresolution microscopy types (for example Structured illumination microscopy (SIM), Stochastic optical reconstruction Microscopy or Stimulated emission depletion (STORM) microscopy (STED)), ExM provides the solution to achieve nano-scale resolution (~70 nm) through the physical expansion of biological samples. independent of an intricate imaging system¹. The fundamental of the physical expansion of biological specimens is to infuse fixed samples with the gelling solution consisting of monomers, followed by in situ polymerization for the formation of a tissue hydrogel hybrid. The gelling solution consists of sodium acrylate, a monomer that confers superabsorbent properties, along with co-monomer acrylamide, cross-linker N-N'-methylenebisacrylamide, accelerator tetramethylethylenediamine, and initiator ammonium persulfate, to trigger the free-radical polymerization. After homogenizing the mechanical properties of the sample-gel composite, isotropic expansion can be achieved after iterative washes with pure water. The extent of expansion, which is defined as expansion factor, can be controlled by varying the salt concentration of the solution, where the maximum expansion factor is reached with pure water.



Figure 1. Nanoscale imaging of E. coli cells by expansion microscopy. (A) Image of pre-expansion *E.coli* TOP10 cells. Stain: DAPI, Blue; WGA, Magenta. (B) Image of expanded *E.coli* cells in the same field of view as A. (C) Zoom-in field of view of A as indicated by white dash box. (D) Zoom-in field of view of B as indicated by white dash box. Scale bar (in biological scale): (A-B) 10 μ m; (C-D) 2 μ m. Expansion factor: 2.8 in 1× PBS.

Parameters in the steps of gelation and homogenization and the strategies for presenting target molecules have been continually improving to accommodate a broader range of applications. The first version of ExM requires the unique design of the fluorescent labels that can directly incorporate to the polymer network to survive the homogenization with protease K^2 . The proteinretention ExM (proExM) later managed to preserve antigens in the gel via modification of the amines on proteins with an acrylamide functional group that participate in the free-radical polymerization for anchoring³. Successful expansions have been performed among a wide range of cell lineages, and tissue types⁴, showing the potential of ExM to facilitate the basic science research and clinical sample diagnosis.

Owing to the small size, microbial cells, such as bacteria and fungi, are challenging to study using conventional optical microscope. Nanoscale optical imaging of microbial cells can potentially open new doors to experimental observation of novel biological structure with molecular contrast, which is not available with electron microscopy. Escherichia coli (*E. coli*) is one of the most studied gram-negative bacteria. *E. coli* possesses a cell wall composed of a thin peptidoglycan layer (2-4 nm) and an outer membrane layer. The peptidoglycan polysaccharide chain is composed alternatively from N-acetylglucosamine and N-acetylmuramic acids and is further crosslinked by tetrapeptides, exhibiting rigidity for maintaining cell shape and providing protection from osmotic lysis⁵. In a recent study, application of ExM on imaging *E. coli* cells has been reported in combination with a miniaturized microscope⁶, using the first generation ExM^2 .

Here, we demonstrate the application of ExM on imaging E. coli cells with ~80 nm resolution. E. coli TOP10 cells, a common chemical competent E. coli strain, were fixed with 4% paraformaldehyde (PFA) and stained with DAPI (4',6-diamidino-2phenylindole) and a lectin wheat germ agglutinin (WGA) to visualize its nucleus and polysaccharidecontaining cell structure, respectively (Figure 1). Images were all acquired with a Nikon Ti2 eclipse fluorescence microscope, with a CSU-W1 confocal module under a long working distance water immersion 40× objective (numerical aperture (NA) = 1.33, confocal microscope). Without expansion it is hard to clearly distinguish single cells and both DAPI and WGA stains overlapped within a rod shape (Figure 1A, 1C). The treatment of enzymatic homogenization allows isotropic expansion of the gelled E. coli cells. Approximately 3-fold expansion factor (2.8) was achieved in PBS (Figure 1B). After expansion, the cell morphology during binary fission is clearly resolved while the overall distortion compared to pre-expansion state is minimal (Figure 1D). Further expansion is achievable with pure water, if additional resolving power is needed.

Acknowledgements

This study was supported by the faculty start up fund from Carnegie Mellon University and NIH Director's New Innovator Award (DP2 OD025926-01) to YZ. We thank Saigopalakrishna Yerneni for providing E. coli samples.

Conflict of Interest

YZ is an inventor in multiple ExM-related patents/ patent applications.

References

- 1. Wassie AT, Zhao Y, Boyden ES. Expansion microscopy: principles and uses in biological research. Nat Methods. 2019; 16(1):33-41. doi:10.1038/s41592-018-0219-4
- Chen F, Tillberg PW, Boyden ES. Expansion Microcopy. Science.2015 Jan 30; 347(6221):543-8. doi:10.1126/science.1260088
- 3. Tillberg PW, Chen F, Piatkevich KD, et al. Protein-retention expansion microscopy of cells and tissues labeled using standard fluorescent proteins and antibodies. Nat Biotechnol. 2016; 34(9):987-992. doi:10.1038/nbt.3625.
- 4. Zhao Y, Bucur O, Irshad H, et al. Nanoscale imaging of clinical specimens using pathologyoptimized expansion microscopy. Nat Biotech. 2017; 35(8):757-764. doi: 10.1038/nbt.3892.
- 5. Costerton JW, Ingram JM, Cheng KJ. Structure and function of the cell envelope of gram negative bacteria. Bacteriol Rev. 1974.
- 6. Zhang YS, Chang J-B, Alvarez MM, et al. Hybrid Microscopy: Enabling Inexpensive High-Performance Imaging through Combined Physical and Optical Magnifications. Sci Rep. 2016;6(1):22691. doi:10.1038/srep22691

DISCOVERIES is a peer-reviewed, open access, online, multidisciplinary and integrative journal, *publishing high impact and innovative manuscripts* from all areas related to MEDICINE, BIOLOGY and CHEMISTRY.

This article is an Open Access article distributed under the terms of the Creative Commons Attribution License (https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited and is not used for commercial purposes; 2019, Applied Systems;